

PRODUCT INFORMATION

Target	CDX1
Synonyms	N/A
Description	Purified recombinant protein of Human caudal type homeobox 1 (CDX1), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Delivery	2-3 weeks
Uniprot ID	P47902
Expression Host	E. coli
Tag	N-His
Molecular Characterization	N/A
Molecular Weight	28 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25mM Tris, pH8.0, 150 mM NaCl, 10% glycerol, 1 % Sarkosyl.
Storage & Shipping	Store at -80°C.
Background	This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded DNA-binding protein regulates intestine-specific gene expression and enterocyte differentiation. It has been shown to induce expression of the intestinal alkaline phosphatase gene, and inhibit beta-catenin/T-cell factor transcriptional activity. [provided by RefSeq, Jul 2008]
Usage	Research use only
Conjugate	Unconjugated

