

## PRODUCT INFORMATION

<b>Target</b>	CLU
<b>Synonyms</b>	AAG4;APOJ;CLI;KUB1;MGC24903;SGP-2;SGP2;SP-40;TRPM-2;TRPM2
<b>Description</b>	Recombinant human CLU protein with C-terminal 6×His tag
<b>Delivery</b>	In Stock
<b>Uniprot ID</b>	P10909
<b>Expression Host</b>	HEK293
<b>Tag</b>	C-6×His Tag
<b>Molecular Characterization</b>	CLU(Asp23-Glu449) 6×His tag
<b>Molecular Weight</b>	The protein has a predicted molecular mass of 50.9 kDa after removal of the signal peptide. The apparent molecular mass of CLU-His is approximately 35-55 kDa due to glycosylation.
<b>Purity</b>	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Storage &amp; Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	The protein encoded by this gene is a secreted chaperone that can under some stress conditions also be found in the cell cytosol. It has been suggested to be involved in several basic biological events such as cell death, tumor progression, and neurodegenerative disorders. Alternate splicing results in both coding and non-coding variants.[provided by RefSeq, May 2011]
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated






Figure 1. Human CLU Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

