

PRODUCT INFORMATION

Target	CNPY2
Synonyms	HP10390; MSAP; TMEM4; ZSIG9
Description	Recombinant protein of human canopy 2 homolog (zebrafish) (CNPY2)
Delivery	2-3 weeks
Uniprot ID	Q9Y2B0
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	20.5 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Positive regulator of neurite outgrowth by stabilizing myosin regulatory light chain (MRLC). It prevents MIR-mediated MRLC ubiquitination and its subsequent proteasomal degradation.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

