

PRODUCT INFORMATION

Target	CRBN
Synonyms	N/A
Description	Purified recombinant protein of Human cereblon (CRBN), transcript variant 1, full length, with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Delivery	2-3 weeks
Uniprot ID	Q96SW2
Expression Host	HEK293T
Tag	Myc-DDK
Molecular Characterization	N/A
Molecular Weight	50.5 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	A DNA sequence from TrueORF clone, [RC209794], encoding the full-length of CRBN
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

