

PRODUCT INFORMATION

Target	DEPTOR
Synonyms	DEP.6; DEPDC6
Description	Recombinant protein of human DEP domain containing 6 (DEPDC6)
Delivery	2-3 weeks
Uniprot ID	Q8TB45
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	46.1 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Negative regulator of the mTORC1 and mTORC2 signaling pathways. Inhibits the kinase activity of both complexes.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMA BIO CONFIDENTIAL

