

PRODUCT INFORMATION

Target	DLL4
Synonyms	AOS6; delta4; hdelta2
Description	Recombinant human DLL4 Protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	Q9NR61
Expression Host	HEK293
Tag	C-6×His tag
Molecular Characterization	DLL4(Ser27-Pro524) 6×His tag
Molecular Weight	The protein has a predicted molecular mass of 55.1 kDa after removal of the signal peptide. The apparent molecular mass of DLL4-His is approximately 55-70 kDa due to glycosylation.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	This gene is a homolog of the Drosophila delta gene. The delta gene family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. [provided by RefSeq, Jul 2008]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human DLL4 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

