

PRODUCT INFORMATION

Target	DNAAF6
Synonyms	CILD36; CXorf41; NYSAR97; PIH1D3; TWISTER
Description	Recombinant protein of human chromosome X open reading frame 41 (CXorf41)
Delivery	2-3 weeks
Uniprot ID	Q9NQM4
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	23.9 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Plays a role in cytoplasmic pre-assembly of axonemal dynein.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

