

PRODUCT INFORMATION

Target	EFCAB4A
Synonyms	EFCAB4A
Description	Recombinant protein of human EF-hand calcium binding domain 4A (EFCAB4A)
Delivery	2-3 weeks
Uniprot ID	Q8N4Y2
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	32.6 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Plays a role in store-operated Ca(2+) entry (SOCE).[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

