

PRODUCT INFORMATION

Target	EMCN
Synonyms	Endomucin;MUC-14
Description	Recombinant human EMCN protein with C-terminal human Fc tag
Delivery	In Stock
Uniprot ID	Q9ULC0
Expression Host	HEK293
Tag	C-Human Fc Tag
Molecular Characterization	EMCN(Asn19-Ser190) hFc(Glu99-Ala330)
Molecular Weight	The protein has a predicted molecular mass of 43.9 kDa after removal of the signal peptide. The apparent molecular mass of EMCN-hFc is approximately 55-100 kDa due to glycosylation.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	EMCN is a mucin-like sialoglycoprotein that interferes with the assembly of focal adhesion complexes and inhibits interaction between cells and the extracellular matrix (Kinoshita et al., 2001 [PubMed 11418125]).[supplied by OMIM, Mar 2008]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human EMCN Protein, hFc Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

