

**PRODUCT INFORMATION**

|   |   |
|---|---|
| <b>Target</b>                           | ENPP3   |
| <b>Synonyms</b>                         | B10; NPP3; PDNP3; CD203c; PD-IBETA  |
| <b>Description</b>                      | Recombinant human ENPP3(546-875) Protein with N-terminal 10×His tag   |
| <b>Delivery</b>                         | In Stock  |
| <b>Uniprot ID</b>                       | O14638  |
| <b>Expression Host</b>                  | HEK293  |
| <b>Tag</b>                              | N-10×His tag  |
| <b>Molecular Characterization</b>       | 10×His tag ENPP3(Pro546-Ile875)   |
| <b>Molecular Weight</b>                 | The protein has a predicted molecular mass of 39.0 kDa after removal of the signal peptide. The apparent molecular mass of His-ENPP3(546-875) is approximately 55-70 kDa due to glycosylation.  |
| <b>Purity</b>                           | The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.  |
| <b>Formulation &amp; Reconstitution</b> | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.   |
| <b>Yefei_Storage</b>                    | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.   |
| <b>Background</b>                       | The protein encoded by this gene belongs to a series of ectoenzymes that are involved in hydrolysis of extracellular nucleotides. These ectoenzymes possess ATPase and ATP pyrophosphatase activities and are type II transmembrane proteins. Expression of the related rat mRNA has been found in a subset of immature glial cells and in the alimentary tract. The corresponding rat protein has been detected in the pancreas, small intestine, colon, and liver. The human mRNA is expressed in glioma cells, prostate, and uterus. Expression of the human protein has been detected in uterus, basophils, and mast cells. Two transcript variants, one protein coding and the other non-protein coding, have been found for this gene. [provided by RefSeq, Oct 2015] |
| <b>Usage</b>                            | Research use only   |
| <b>Conjugate</b>                        | Unconjugated  |






Figure 1. Human ENPP3(546-875) Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

