

**PRODUCT INFORMATION**

<b>Target</b>	FAM170A
<b>Synonyms</b>	ZNFD
<b>Description</b>	Purified recombinant protein of Human family with sequence similarity 170, member A (FAM170A), transcript variant 1, full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
<b>Delivery</b>	2-3 weeks
<b>Uniprot ID</b>	A1A519
<b>Expression Host</b>	E. coli
<b>Tag</b>	N-His
<b>Molecular Characterization</b>	N/A
<b>Molecular Weight</b>	36.9 kDa
<b>Purity</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Formulation &amp; Reconstitution</b>	25mM Tris, pH8.0, 150 mM NaCl, 10% glycerol, 1 % Sarkosyl
<b>Storage &amp; Shipping</b>	Store at -80°C.
<b>Background</b>	Acts as a nuclear transcription factor that positively regulates the expression of heat shock genes. Binds to heat shock promoter elements (HSE).[UniProtKB/Swiss-Prot Function]
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

