

PRODUCT INFORMATION

Target	FOXN2
Synonyms	HTLF
Description	Recombinant protein of human forkhead box N2 (FOXN2)
Delivery	2-3 weeks
Uniprot ID	P32314
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	47 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	This gene encodes a forkhead domain binding protein and may function in the transcriptional regulation of the human T-cell leukemia virus long terminal repeat. [provided by RefSeq, Jul 2008]
Usage	Research use only
Conjugate	Unconjugated

