

**PRODUCT INFORMATION**

<b>Target</b>	FOXR2
<b>Synonyms</b>	FOXN6
<b>Description</b>	Purified recombinant protein of Human forkhead box R2 (FOXR2), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
<b>Delivery</b>	1 week
<b>Uniprot ID</b>	Q6PJQ5
<b>Expression Host</b>	E. coli
<b>Tag</b>	N-His
<b>Molecular Characterization</b>	N/A
<b>Molecular Weight</b>	35.7 kDa
<b>Purity</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Formulation &amp; Reconstitution</b>	50mM Tris, 8M Urea, pH8.0.
<b>Storage &amp; Shipping</b>	Store at -80°C.
<b>Background</b>	A DNA sequence encoding human full-length FOXR2
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

DIMABIO CONFIDENTIAL

