

PRODUCT INFORMATION

Target	GABPB2
Synonyms	GABPB-2
Description	Recombinant protein of human GA binding protein transcription factor, beta subunit 2 (GABPB2)
Delivery	2-3 weeks
Uniprot ID	Q8TAK5
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	48.5 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	May function as transcription factor capable of interacting with purine rich repeats (GA repeats).[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

