

PRODUCT INFORMATION

Target	GBP4
Synonyms	Mpa2
Description	Recombinant protein of human guanylate binding protein 4 (GBP4)
Delivery	2-3 weeks
Uniprot ID	Q96PP9
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	73 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Guanylate-binding proteins, such as GBP4, are induced by interferon and hydrolyze GTP to both GDP and GMP (Vestal, 2005 [PubMed 16108726]).[supplied by OMIM, Dec 2008]
Usage	Research use only
Conjugate	Unconjugated

