

PRODUCT INFORMATION

Tag	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 C-Flag Tag
Target	GBRR3
Synonyms	N/A
Description	Human GBRR3 full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	A8MPY1
Expression Host	HEK293
Protein Families	Ion Channels: Cys-loop Receptors
Protein Pathways	N/A
Molecular Weight	The human full length GBRR3 protein has a MW of 54.3kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	The neurotransmitter gamma-aminobutyric acid (GABA) functions in the central nervous system to regulate synaptic transmission of neurons. This gene encodes one of three related subunits, which combine as homo- or hetero-pentamers to form GABA(C) receptors. In humans, some individuals contain a single-base polymorphism (dbSNP rs832032) that is predicted to inactivate the gene product. [provided by RefSeq, Jan 2012]
Usage	Research use only
Conjugate	Unconjugated

