

PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806 C-Hag Tag Tag

Synonyme

Human GPR45 full length protein-synthetic nanodisc

Description Delivery 6~8weeks Uniprot ID O9Y5Y3 HEK293 Expression Host

Protein Families GPCR, Transmembrane, Druggable Genome,

GPCRDB Class A Rhodopsin-like, Protein Pathways Molecular Weight

The human full length GPR45 protein has a MW of 42kDa

DIMARIO CONTINUE DE LA CONTINUE DE L Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCI, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

ryup-mizauum, r-rease see Letrincate or Analysis for Specific Instructions of reconstitution. Store at -20°C to -80°C for 12 months in Jophilized form. After reconstitution, in fon intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. This intronless gene encodes a member of the G protein-coupled receptor (GPCR) family. Members of this protein family contain seven putative transmembrane domains and may mediate signaling processes to the interior of the cell via activation of heterotrimeric G proteins. This protein may function in the central nervous system, [provided by keFSeq. Jul 2008] Storage & Shipping

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only Conjugate

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

