

PRODUCT INFORMATION

Warning. Undefined variable ShasAttributeValueDescription in Chwwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 CFlagSttrp En Tag

Synonyme

Human GPRC5D-Strep full length protein-synthetic nanodisc

Description Delivery 6~8weeks Uniprot ID O9NZD1 HEK293 Expression Host Protein Families GPCR Protein Pathways

Molecular Weight The human full length GPRC5D-Strep Protein has a MW of 38.6 kDa

Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCI, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

syspimication: recease see Lettinicate or Analysis for Specific Instructions of reconstitution. Store at -20°C to -80°C for 12 months in Joyahitzed form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

The protein encoded by this gene is a member of the G protein-coupled receptor family. Recent studies demonstrate that GPRC5D is expressed on malignant bone marrow plasma cells, whereas normal tissue expression is limited to the hair follicle. It may represent a potential target for effector-cell-mediated therapy to treat plasma-cell disorders like MM. Storage & Shipping

Email: info@dimabio.com Website: www.dimabio.com

Usage Conjugate

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

