

PRODUCT INFORMATION

Tag	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 C-Flag&Strep Tag
Target	GRL1
Synonyms	GCOM1, GRINL1A, Gdown, Gdown1
Description	Human GRL1-Strep full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	POCAP2
Expression Host	HEK293
Protein Families	Ion Channels: Other
Protein Pathways	N/A
Molecular Weight	The human full length GRL1-Strep protein has a MW of 41.7 kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	This gene encodes a subunit of a specific form of RNA polymerase II termed Pol II(G). The encoded protein may act as a negative regulator of transcriptional activation by the Mediator complex. Alternative splicing results in multiple transcript variants. There is a pseudogene for this gene on chromosome 4. Readthrough transcription between this gene and the neighboring upstream gene MYZAP (myocardial zonula adherens protein) is represented with GenelD 145791. (provided by RefSeq, Oct 2013)
Usage	Research use only
Conjugate	Unconjugated

