

PRODUCT INFORMATION

Target	GTPBP10
Synonyms	ObgH2;µg0751c10
Description	Recombinant protein of human GTP-binding protein 10 (putative) (GTPBP10), transcript variant 2
Delivery	2-3 weeks
Uniprot ID	A4D1E9
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	42.8 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Small G proteins, such as GTPBP10, act as molecular switches that play crucial roles in the regulation of fundamental cellular processes such as protein synthesis, nuclear transport, membrane trafficking, and signal transduction (Hirano et al., 2006 [PubMed 17054726]).[supplied by OMIM, Mar 2008]
Usage	Research use only
Conjugate	Unconjugated

