

**PRODUCT INFORMATION**

<b>Target</b>	HDDC2
<b>Synonyms</b>	C6orf74; CGI-130; dj16705.2; NS5ATP2
<b>Description</b>	Recombinant protein of human HD domain containing 2 (HDDC2)
<b>Delivery</b>	2-3 weeks
<b>Uniprot ID</b>	Q7Z4H3
<b>Expression Host</b>	HEK293T
<b>Tag</b>	C-Myc/DDK
<b>Molecular Characterization</b>	N/A
<b>Molecular Weight</b>	23.2 kDa
<b>Purity</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Formulation &amp; Reconstitution</b>	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
<b>Storage &amp; Shipping</b>	Store at -80°C.
<b>Background</b>	Catalyzes the dephosphorylation of the nucleoside 5'-monophosphates deoxyadenosine monophosphate (dAMP), deoxycytidine monophosphate (dCMP), deoxyguanosine monophosphate (dGMP) and deoxythymidine monophosphate (dTMP).[UniProtKB/Swiss-Prot Function]
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

