

PRODUCT INFORMATION

Target	HSDL2
Synonyms	C9orf99; SDR13C1
Description	Recombinant protein of human hydroxysteroid dehydrogenase like 2 (HSDL2)
Delivery	2-3 weeks
Uniprot ID	Q6YN16
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	45.2 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Has apparently no steroid dehydrogenase activity.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

