

PRODUCT INFORMATION

IL19 Target

Synonyms IL-10C; MDA1; NG.1; ZMDA1

Recombinant Human IL19 Protein with C-terminal **Description**

human Fc tag

Delivery In Stock **Uniprot ID** Q9UHD0 **Expression Host HEK293**

Tag C-Human Fc Tag

Molecular

Background

IL19(Leu25-Ala177) hFc(Glu99-Ala330) Characterization

The protein has a predicted molecular mass of **Molecular Weight**

43.9 kDa after removal of the signal peptide. The apparent molecular mass of IL19-hFc is

approximately 35-70 kDa due to glycosylation. The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue

Purity

staining.

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution. Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

temperature.

The protein encoded by this gene is a cytokine that belongs to the IL10 cytokine subfamily. This cytokine is found to be preferentially expressed in monocytes. It can bind the IL20 receptor complex and lead to the activation of the signal transducer and activator of transcription 3 (STAT3). A similar cytokine in mouse is reported to up-regulate the expression of IL6 and TNF-alpha and induce


apoptosis, which suggests a role of this cytokine in inflammatory responses. Alternatively spliced transcript variants encoding the distinct isoforms have been described. [provided by RefSeq, Jul

2008]

Usage Research use only

Conjugate Unconjugated













Figure 1.Human IL19 Protein, hFc Tag on SDS-PAGE under reducing condition.

Human IL19, hFc Tagged protein ELISA

0.2 μg of Human IL19, hFc tagged protein per well



Biotinylated Anti-IL19 antibody, Rabbit mAb (ng/mL)

Figure 2. ELISA plate pre-coated by 2 μ g/mL (100 μ L/well) Human IL19 Protein, hFc Tag (PME101250) can bind Biotinylated Anti-IL19 antibody, Rabbit mAb (DME101184B) in a linear range of 16–80 ng/mL.

Email: info@dimabio.com Website: www.dimabio.com

