

PRODUCT INFORMATION

Target	IRGC
Synonyms	CINEMA; IFGGE; ligp5; IRGC1; R30953_1
Description	Recombinant protein of human immunity-related GTPase family, cinema (IRGC)
Delivery	2-3 weeks
Uniprot ID	Q6NXR0
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	50.1 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris-HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Recombinant protein was produced with TrueORF clone, [RC209793]. Click on the TrueORF clone link to view cDNA and protein sequences.
Usage	Research use only
Conjugate	Unconjugated

