

PRODUCT INFORMATION

Target	ITGB1
Synonyms	GPIIA;CD29;VLA-4 subunit beta
Description	Recombinant human ITGB1 protein with C-terminal human Fc tag
Delivery	In Stock
Uniprot ID	P05556
Expression Host	HEK293
Tag	C-Human Fc Tag
Molecular Characterization	ITGB1(Gln21-Asp728) hFc(Glu99-Ala330)
Molecular Weight	The protein has a predicted molecular mass of 104.5 kDa after removal of the signal peptide. The apparent molecular mass of ITGB1-hFc is approximately 130-250 kDa due to glycosylation.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Integrins are heterodimeric proteins made up of alpha and beta subunits. At least 18 alpha and 8 beta subunits have been described in mammals. Integrin family members are membrane receptors involved in cell adhesion and recognition in a variety of processes including embryogenesis, hemostasis, tissue repair, immune response and metastatic diffusion of tumor cells. This gene encodes a beta subunit. Multiple alternatively spliced transcript variants which encode different protein isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human IL17RB Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

