

PRODUCT INFORMATION

Tag	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 C-Flag&Strep Tag
Target	KCNE5
Synonyms	KCNE1L
Description	Human KCNE5-Strep full length protein-synthetic nanodisc
Delivery	6-8weeks
Uniprot ID	Q9UJ90
Expression Host	HEK293
Protein Families	Ion Channels: Other
Protein Pathways	N/A
Molecular Weight	The human full length KCNE5-Strep protein has a MW of 15 kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	This gene encodes a member of a family of single pass transmembrane domain proteins that function as ancillary subunits to voltage-gated potassium channels. Members of this family affect diverse processes in potassium channel regulation, including ion selectivity, voltage dependence, and anterograde recycling from the plasma membrane. Variants of this gene are associated with idiopathic ventricular fibrillation and Brugada syndrome. [provided by RefSeq, Nov 2016]
Usage	Research use only
Conjugate	Unconjugated

