

PRODUCT INFORMATION

Target	KIAA0247
Synonyms	DRAGO; KIAA0247
Description	Recombinant protein of human KIAA0247 (KIAA0247)
Delivery	2-3 weeks
Uniprot ID	Q92537
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	31.9 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	May play a role in growth-suppressive activity and cell death (PubMed:24652652). May be involved in the production of chemokine molecules in umbilical vein endothelial cells (HUVECs) cultured in THP1 monocyte LPS-induced medium (PubMed:20236627). Plays a role in preventing tumor onset (By similarity).[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

