

PRODUCT INFORMATION

Target	KREMEN2
Synonyms	KRM2
Description	Recombinant Human KREMEN2 Protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	Q8NCW0
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	KREMEN2(Gly26-Ala364) 6×His tag
Molecular Weight	The protein has a predicted molecular mass of 36.8 kDa after removal of the signal peptide.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	This gene encodes a high-affinity dickkopf homolog 1 (DKK1) transmembrane receptor. A similar protein in mouse functions interacts with DKK1 to block wingless (WNT)/beta-catenin signaling. The encoded protein forms a ternary membrane complex with DKK1 and the WNT receptor lipoprotein receptor-related protein 6 (LRP6), and induces rapid endocytosis and removal of LRP6 from the plasma membrane. It contains extracellular kringle, WSC, and CUB domains. Alternatively spliced transcript variants encoding distinct isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human KREMEN2 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

