

PRODUCT INFORMATION

Target	LCN6
Synonyms	hLcn5; LCN5; UNQ643
Description	Recombinant protein of human lipocalin 6 (LCN6)
Delivery	1 week
Uniprot ID	P62502
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	17.9 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	May play a role in male fertility.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

