

PRODUCT INFORMATION

Target	MAN2B2
Synonyms	N/A
Description	Recombinant protein of human mannosidase, alpha, class 2B, member 2 (MAN2B2)
Delivery	2-3 weeks
Uniprot ID	Q9Y2E5
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	113.8 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Recombinant protein was produced with TrueORF clone, [RC210523]. Click on the TrueORF clone link to view cDNA and protein sequences.
Usage	Research use only
Conjugate	Unconjugated

DIMA BIO CONFIDENTIAL

