

PRODUCT INFORMATION

Target	METTL4
Synonyms	HsT661
Description	Recombinant protein of human methyltransferase like 4 (METTL4)
Delivery	2-3 weeks
Uniprot ID	Q8N3J2
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	53.9 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Probable methyltransferase.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

