

PRODUCT INFORMATION

Target	MFAP1
Synonyms	AMP
Description	Purified recombinant protein of Human microfibrillar-associated protein 1 (MFAP1), with N-terminal HIS tag, expressed in E.Coli, 50ug
Delivery	2-3 weeks
Uniprot ID	P55081
Expression Host	E. coli
Tag	N-His
Molecular Characterization	N/A
Molecular Weight	51.8 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25mM Tris, pH8.0, 150 mM NaCl, 10% glycerol, 1 % Sarkosyl.
Storage & Shipping	Store at -80°C.
Background	Involved in pre-mRNA splicing as a component of the spliceosome.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

