

PRODUCT INFORMATION

Target	MIOX
Synonyms	ALDRL6
Description	Recombinant protein of human myo-inositol oxygenase (MIOX)
Delivery	2-3 weeks
Uniprot ID	Q9UGB7
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	32.8 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Recombinant protein was produced with TrueORF clone, [RC210070]. Click on the TrueORF clone link to view cDNA and protein sequences.
Usage	Research use only
Conjugate	Unconjugated

DIMA BIO CONFIDENTIAL

