

PRODUCT INFORMATION

Target	BMP-2
Synonyms	Bone Morphogenetic Protein 2;BMP-2;Bone Morphogenetic Protein 2A;BMP-2A;BMP2;BMP2A
Description	Recombinant Human/Mouse/Rat Bone Morphogenetic Protein 2 is produced by our E.coli expression system and the target gene encoding Gln283-Arg396 is expressed.
Delivery	In Stock
Uniprot ID	P12643
Expression Host	E.coli
Tag	
Molecular Characterization	Not available
Molecular Weight	13.3 KDa
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Formulation & Reconstitution	Lyophilized from a 0.2 μ m filtered solution of 10mM HAC-NH ₄ Ac, 4% D-Mannitol, pH 4.0.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Bone Morphogenetic Protein-2 (BMP-2) is one of the bone-growth regulatory factors that belong to the transforming growth factor-beta (TGF-beta) superfamily of proteins. BMPs are synthesized as large precursor molecules, which are cleaved by proteolytic enzymes. The active form of BMP-2 can consist of a dimer of two identical proteins or a heterodimer of two related bone morphogenetic proteins.
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Greater than 95% as determined by reducing SDS-PAGE.

DIMABIO CONFIDENTIAL

