

PRODUCT INFORMATION

Target	MSP1D1
Synonyms	APOA1
Description	Recombinant human MSP1D1(del H5) Protein with N-terminal 6×His tag
Delivery	In Stock
Uniprot ID	P02647
Expression Host	HEK293
Tag	N-6×His tag
Molecular Characterization	6×His tag APOA1(Ser79-Gln267) del(Pro145-Ser166)
Molecular Weight	The protein has a predicted molecular mass of 21.2 kDa after removal of the signal peptide. The apparent molecular mass of His-MSP1D1(del H5) is approximately 15-25 kDa due to glycosylation.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	MSP1D1 (Membrane Scaffold Protein 1D1) is derived from Apolipoprotein A-I (ApoA-I) and is used in the formation of nanodiscs, which are tools for studying membrane proteins in a controlled lipid environment. The "H" in H5 refers to a histidine residue at position 5 in the protein sequence. The "del" symbol indicates a deletion.
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human MSP1D1(del H5) Protein, His Tag on SDS-PAGE under reducing condition.

