

PRODUCT INFORMATION

Target	N4BP3
Synonyms	LZTS4
Description	Recombinant protein of human Nedd4 binding protein 3 (N4BP3)
Delivery	1 week
Uniprot ID	O15049
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	60.3 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Plays a role in axon and dendrite arborization during cranial nerve development. May also be important for neural crest migration and early development of other anterior structures including eye, brain and cranial cartilage.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

