

PRODUCT INFORMATION

Target	N/A
Synonyms	N/A
Description	Purified recombinant protein of Human nudix (nucleoside diphosphate linked moiety X)-type motif 8 (NUDT8), full length, with N-terminal HIS tag, expressed in E. coli, 50ug
Delivery	2-3 weeks
Uniprot ID	N/A
Expression Host	E. coli
Tag	N-His
Molecular Characterization	N/A
Molecular Weight	14.9 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25mM Tris, pH8.0, 150mM NaCl, 10% glycerol, 1% Sarkosyl.
Storage & Shipping	Store at -80°C.
Background	A DNA sequence encoding human full-length NUDT8
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

