

**PRODUCT INFORMATION**

<b>Target</b>	NUP42
<b>Synonyms</b>	CG1; hCG1; NLP-1; NLP_1; NUPL2
<b>Description</b>	Purified recombinant protein of Human nucleoporin like 2 (NUPL2), full length, with N-terminal HIS tag, expressed in E. coli, 50ug
<b>Delivery</b>	2-3 weeks
<b>Uniprot ID</b>	O15504
<b>Expression Host</b>	E. coli
<b>Tag</b>	N-His
<b>Molecular Characterization</b>	N/A
<b>Molecular Weight</b>	44.7 kDa
<b>Purity</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Formulation &amp; Reconstitution</b>	50mM Tris, 8M Urea, pH8.0
<b>Storage &amp; Shipping</b>	Store at -80°C. Required for the export of mRNAs containing poly(A) tails from the nucleus into the cytoplasm. In case of infection by HIV-1, it may participate in the docking of viral Vpr at the nuclear envelope.[UniProtKB/Swiss-Prot Function]
<b>Background</b>	
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

