

PRODUCT INFORMATION

Target	OLIG3
Synonyms	Bhlhb7; bHLHe20
Description	Recombinant protein of human oligodendrocyte transcription factor 3 (OLIG3)
Delivery	2-3 weeks
Uniprot ID	Q7RTU3
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	29.2 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

