

PRODUCT INFORMATION

PDL1 **Target**

PD-L1;CD274;B7-H1;PDCD1L1;PDCD1LG1 **Synonyms** Recombinant human PD-L1 protein with C-

Description terminal mouse Fc and 6×His tag

Delivery In Stock **Uniprot ID** Q9NZQ7 **Expression Host HEK293**

Tag C-Mouse Fc and 6×His Tag

Molecular PD-L1(Phe19-Arg238) mFc(Pro99-Lys330) 6×His

Characterization

Purity

Background

Molecular Weight 52.3 kDa after removal of the signal peptide.

The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue

The protein has a predicted molecular mass of

staining.

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before Formulation & Reconstitution

lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Storage & Shipping

Lyophilized proteins are shipped at ambient

temperature.

This gene encodes an immune inhibitory receptor ligand that is expressed by hematopoietic and non-hematopoietic cells, such as T cells and B cells and various types of tumor cells. The encoded protein is a type I transmembrane protein that has immunoglobulin V-like and C-like domains. Interaction of this ligand with its receptor inhibits T-cell activation and cytokine production. During infection or inflammation of

normal tissue, this interaction is important for

preventing autoimmunity by maintaining homeostasis of the immune response. In tumor microenvironments, this interaction provides an immune escape for tumor cells through cytotoxic T-cell inactivation. Expression of this gene in tumor cells is considered to be prognostic in many types of human malignancies, including

colon cancer and renal cell carcinoma. Alternative splicing results in multiple transcript variants.

Usage Research use only Conjugate Unconjugated

> Email: info@dimabio.com Website: www.dimabio.com










Figure 1. Human PD-L1 Protein, mFc-His Tag on SDS-PAGE under reducing condition.

Human PDL1, mFc-His Tagged protein ELISA

0.2 µg of Human PDL1, mFc-His Tagged protein per well




Figure 2. ELISA plate pre-coated by 2 μ g/ml (100 μ l/well) Human PD-L1, mFc-His tagged protein PME100023 can bind Anti-PDL1 Neutralizing antibody BME100009 in a linear range of 0.24-7.81 ng/ml.

Email: info@dimabio.com Website: www.dimabio.com



Human PD-L1, mFc-His Tagged protein ELISA

0.2 µg of PD1, hFc-His Tagged protein per well




Figure 3. ELISA plate pre-coated by 2 μ g/ml (100 μ l/well) Human PD1, hFc-His tagged protein PME100462 can bind Human PDL1, mFc-His tagged protein (PME100023) in a linear range of 62.5-251.1 ng/ml.




Figure 4. Flow cytometry analysis with 15 μ g/mL Human PD-L1 Protein, mFc-His tag (PME100023) on HEK293 cells transfected with human PD1 (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

Email: info@dimabio.com Website: www.dimabio.com

