

PRODUCT INFORMATION

Target	PRKRIR
Synonyms	DAP4; P52rIPK; PRKRIR; THAP0
Description	Recombinant protein of human protein-kinase, interferon-inducible double stranded RNA dependent inhibitor, repressor of (P58 repressor) (PRKRIR)
Delivery	2-3 weeks
Uniprot ID	O43422
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	87.5 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Upstream regulator of interferon-induced serine/threonine protein kinase R (PKR). May block the PKR-inhibitory function of DNAJC3, resulting in restoration of kinase activity and suppression of cell growth.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

