

PRODUCT INFORMATION

Target	PRTN3
Synonyms	ACPA;AGP7;C-ANCA;CANCA;MBN;MBT;NP-4;NP4;P29;PR-3;PR3
Description	Recombinant Human PRTN3 Protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	P24158
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	PRTN3(Ala26-Pro256) 6×His tag
Molecular Weight	The protein has a predicted molecular mass of 26.2 kDa after removal of the signal peptide. The apparent molecular mass of PRTN3-His is approximately 25-35 kDa due to glycosylation.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Serine protease that degrades elastin, fibronectin, laminin, vitronectin, and collagen types I, III, and IV (in vitro) (PubMed:3198760, PubMed:2033050, PubMed:28240246). By cleaving and activating receptor F2RL1/PAR-2, enhances endothelial cell barrier function and thus vascular integrity during neutrophil transendothelial migration (PubMed:23202369). May play a role in neutrophil transendothelial migration, probably when associated with CD177 (PubMed:22266279).[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human PRTN3 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

