

PRODUCT INFORMATION

Target	RAB26
Synonyms	V46133
Description	Recombinant protein of human RAB26, member RAS oncogene family (RAB26)
Delivery	2-3 weeks
Uniprot ID	Q9ULW5
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	27.7 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Members of the RAB protein family, including RAB26, are important regulators of vesicular fusion and trafficking. The RAB family of small G proteins regulates intercellular vesicle trafficking, including exocytosis, endocytosis, and recycling (summary by Seki et al., 2000 [PubMed 11043516]).[supplied by OMIM, Nov 2010]
Usage	Research use only
Conjugate	Unconjugated

