

PRODUCT INFORMATION

Target	RAB6C
Synonyms	WTH3
Description	Purified recombinant protein of Human RAB6C, member RAS oncogene family (RAB6C), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Delivery	2-3 weeks
Uniprot ID	Q9H0N0
Expression Host	E. coli
Tag	N-His
Molecular Characterization	N/A
Molecular Weight	28.2 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	50mM Tris,8M Urea,pH8.0.
Storage & Shipping	Store at -80°C.
Background	May be involved in the regulation of centrosome duplication and cell cycle progression.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

