

PRODUCT INFORMATION

Target	Rab9
Synonyms	RAB9
Description	Recombinant protein of human RAB9A, member RAS oncogene family (RAB9A)
Delivery	2-3 weeks
Uniprot ID	P51151
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	22.7 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Involved in the transport of proteins between the endosomes and the trans Golgi network. Involved in the recruitment of SGSM2 to melanosomes and is required for the proper trafficking of melanogenic enzymes TYR, TYRP1 and DCT/TYRP2 to melanosomes in melanocytes.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

