

PRODUCT INFORMATION

Target	RANBP6
Synonyms	N/A
Description	Recombinant protein of human RAN binding protein 6 (RANBP6), transcript variant 1
Delivery	2-3 weeks
Uniprot ID	O60518
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	124.5 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	May function in nuclear protein import as nuclear transport receptor.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

