

PRODUCT INFORMATION

Target	RHOJ
Synonyms	ARHJ; RASL7B; TC10B; TCL
Description	Recombinant protein of human ras homolog gene family, member J (RHOJ)
Delivery	2-3 weeks
Uniprot ID	Q9H4E5
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	23.6 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	This gene encodes one of the many small GTP-binding proteins in the Rho family shown to be associated with focal adhesions in endothelial cells (PMID: 21148427, 22103495). The encoded protein is activated by vascular endothelial growth factor and may regulate angiogenesis. [provided by RefSeq, Dec 2011]
Usage	Research use only
Conjugate	Unconjugated

