

PRODUCT INFORMATION

Target	RPL7L1
Synonyms	dj475N16.4
Description	Recombinant protein of human ribosomal protein L7-like 1 (RPL7L1)
Delivery	2-3 weeks
Uniprot ID	Q6DKI1
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	28.5 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Recombinant protein was produced with TrueORF clone, [RC210095]. Click on the TrueORF clone link to view cDNA and protein sequences.
Usage	Research use only
Conjugate	Unconjugated

