

PRODUCT INFORMATION

Target SEZ6 **BSRPC Synonyms**

Recombinant human SEZ6(833-925) Protein with **Description**

C-terminal human Fc tag

Delivery In Stock **Uniprot ID** Q53EL9 **Expression Host HEK293**

Tag C-Human Fc tag

Molecular

Purity

SEZ6(Glu833-His925) hFc(Glu99-Ala330) Characterization

The protein has a predicted molecular mass of **Molecular Weight** 35.8 kDa after removal of the signal peptide.

The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue

staining.

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before Formulation &

lyophilization. Please see Certificate of Analysis Reconstitution

for specific instructions of reconstitution. Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Storage & Shipping

Lyophilized proteins are shipped at ambient

temperature.

The protein encoded by this gene is thought to contain five cysteine-rich motifs that are similar to sushi domains, as well as two domains similar to the amino terminal half of the CUB (for

Background

complement C1r/C1s, Uegf, Bmp1) domain. Mutations in this gene have been associated with febrile seizures. [provided by RefSeq, Jul 2016]

> Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only

Conjugate Unconjugated






Figure 1. Human SEZ6(833-925) Protein, hFc Tag on SDS-PAGE under reducing condition.



