

PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806 C-flagsStrep Tag Tag

Synonyme

Human SMO-Strep full length protein-synthetic nanodisc

Description Delivery 6~8weeks Uniprot ID 099835 HEK293 Expression Host

Protein Families ES Cell Differentiation/IPS ,Transmembrane,Stem cell relevant signaling - DSL/Notch pathway,Druggable Genome,

GPCRDB Other,Hedgehog Netpath 10,Hedgehog Netpath 10,Apoptosis,Cancer,Notch, Protein Pathways

Molecular Weight The human full length SMO-Strep protein has a MW of 86.4 kDa

Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCI, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

DIMARIO CONTINUE DE LA CONTINUE DE L ryupimication: recease see Certificate or Analysis for specific instructions of reconstitution. For the other decisions are a controlled for a see that the controlled for the controlled for a see that the controlled Background

> Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only Conjugate Unconjugated

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

